



# IRON NITRO-PAPS

## Colorimetric determination of iron in serum without blank sample

### TEST SUMMARY

Iron separated from Transferrin, is reduced to bivalent Iron, which reacts with the cromogen present in the reagent, forming a blue compound, which intensity is proportional to the iron concentration present in the sample.

### SAMPLES

Unhemolized serum.  
Stability: 3 days at 2-8°C.

### REAGENTS

Reagent A: Acetate buffer 0.2 mol/l pH 4.3, sodium thioglycolate 40 mmol/l, surfactants.

Reagent B: Nitro-PAPS 0.1 mmol/l pH 4.3.

Standard: lone ferrico 100 µg/dl; stabilizzanti e conservanti.

### MATERIAL REQUIRED BUT NOT SUPPLIED

Normal laboratory equipment. Spectrophotometer UV/VIS with thermostataion. Automatic Micropipette. Cuvette in optical glass or monouse in optical polystyrene. Physiologic solution.

### PRECAUTIONS

Reagent may contain not reactive and conservative components. It is opportune to avoid contacts with the skin and do not swallow. Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

### REAGENTS PREPARATION

Mix a part of Reagent A with a part of Reagent B. Reagents are stable until expiration date on label, stored at 2-8°C.  
Work Reagent is stable 12 days at 2-8°C.

### PROCEDIMENTO

Kind of analysis: End point  
Reading time: 10 minutes  
Wavelength: 582 nm (578-605)  
Temperature: R.T.  
Colour stability: 30 minutes  
Lightpath: 1 cm  
Zero: Blank Reagent

### EXPECTED VALUES

Men	59 - 158 µg/dl
Women	37 - 145 µg/dl

Every laboratory should establish own reference intervals in accordance with own population.

### NOTES

- Particularly turbid sera need the sample blank feasible with only Reagent A.
- Due to the high sensibility of the Reagent, use glassware surely without iron traces.
- If the results are incompatible with clinical presentation, they have to be evaluated within a total clinical study.
- Only for IVD use.

### CALIBRATION/ QUALITY CONTROL

It is suggested to perform an internal quality control using control serum with known iron values.

### TEST PERFORMANCE

#### Precision

Intra-assay (n = 30)	Mean (µg/dl)	SD (µg/dl)	CV%
Sample 1	113.70	1.2905	1.14
Sample 2	165.93	0.8276	0.50

Inter-assay (n = 30)	Mean (µg/dl)	SD (µg/dl)	CV%
Sample 1	113.86	1.3829	1.21
Sample 2	166.60	1.3544	0.81

#### Linearity

The method is linear up to 500 µg/dl

#### Interferences

Sera strong lipemic can interfere with the analysis; it is suggested centrifuge or filtrate the sera with membrane 0.2 µm.

#### Methods comparison

A comparison with a commercial available product gave the following results in a comparison on 30 samples:

Iron Nitro-paps LTA = x  
Iron Nitro-paps competitor = y  
n = 30

$$y = 0,99527x - 0,65775 \quad r = 0,99642$$

#### WASTE DISPOSAL

Product is intended for professional laboratories. Waste products must be handled as per relevant security cards and local regulations.

#### PACKAGING

<b>CODE CC01500 (300 TESTS)</b>	
Reagent A	3 x 100 ml (liquid)
Reagent B	3 x 100 ml (liquid)
Standard	1 x 5 ml (liquid)

### REFERENCES

Weippl.G., et al, Blut. 27, 261 (1973).  
Makino A., Kiyonaga M., Kina K. Clin. Chem. Acta 171:19-28 (1988).  
Maringoni A., Federici G. Eurologo '89 Abstract - Biochimica Clinica. Suppl. 1/8 13,89 178-79.

### MANUFACTURER

LTA s.r.l.  
Via Milano 15/F  
20060 Bussero (Milan) ITALY  
Tel: ++39 02 95409034  
Fax: ++39 02 95334185  
e-mail: info@ltaonline.it  
Website: http://www.ltaonline.it

### SYMBOLS

- Only for IVD use
- Lot of manufacturing
- Code number
- Storage temperature interval
- Expiration date
- Warning, read enclosed documents
- Read the directions
- Biological risk

Mod. 01.06 (ver. 1.4 - 12/02/2009)



REAGENTS	BLANK	STANDARD	SAMPLE
Distilled water	100 µl	--	--
Standard	--	100 µl	--
Sample	--	--	100 µl
Work Reagent	2 ml	2 ml	2 ml

Mix and wait 10 minutes, read the absorbances against Blank at 582 nm.

### CALCULATION

$$\text{Iron } \mu\text{g/dl} = \frac{A_{(\text{sample})}}{A_{(\text{standard})}} \times 100$$