PATHOLOGICAL CONTROL SERUM

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Bibliografia

- 1) Occupational Safety and Health Standards: bloodborne pathogens. (29 CFR 1910.1030). Federal Register. July 1, 1998; 6:267-280
- 2) Council Directive (90/679/EEC). Official Journal of the European Communities. No. L374 from Dec. 31, 1990: 1-12.
- 3) Dati del fabbricante

LTA srl
via Milano 15/F - 20060
Bussero (MI) - ITALY
Tel.++39 02 95409034
Fax.++39 02 95334185

PATHOLOGICAL CONTROL SERUM

UNIVERSAL CONTROL SERUM HUMAN BASED

10 x 5 ml

 CC03200

Intended use

Pathological Control Serum is for use in the quality control of general clinical chemistry methods for the quantitative determination of substrates, electrolytes, lipids, enzymes and proteins. The control is used for monitoring accuracy or precision both for manual techniques and assays on automated clinical chemistry analyzers.

Summary

Pathological Control Serum is a lyophilized control serum based on human serum with concentrations / activities found in the pathological range.

Kit components

Lyophilized control serum 10 x 5 ml

Reactive lyophilizate components

Human serum with chemical additives and tissue extracts of human and animal origin. The origin of biological additives is as follows:

Acid phosphatase	human prostate / potato
Alkaline phosphatase	placenta (human, recombinant)
ALT/GPT	porcine heart
AST/GOT	porcine heart
Albumin	bovine plasma
Aldolase	rabbit muscle
Amylase, total	human saliva / porcine pancreas
Amylase, pancreatic	porcine pancreas
Cholesterol	bovine plasma
Cholinesterase	human serum
Creatine kinase	rabbit muscle
Gamma-GT	porcine kidney
LDH	porcine heart
Lipase	porcine pancreas
Transferrin	human Cohn IV fraction
Triglycerides	chicken egg yolk

The concentrations/activities are lot specific. The exact values and ranges are listed in the enclosed values sheet.

Precautions and warnings

For in vitro diagnostic use only.

Exercise the normal precautions required for handling all laboratory reagents. This product has been prepared exclusively from the blood of donors tested individually and shown by FDA approved methods to be free from HBsAg and antibodies to HCV and HIV. However as no test method can rule out the potential risk of infection with absolute certainty, the material should be treated just as carefully as a patient sample. In the event of exposure the directive of the responsible health authorities should be followed.

Disposal

This product is made to be used in professional laboratories.

Please consult local regulations for a correct waste disposal.

S56: Dispose of this material and its container at hazardous or special waste collection point.

S57: Use appropriate container to avoid environmental contamination.

S61: Avoid release in environment. Refer to special instructions/Safety data sheets.

Handling

Carefully open one bottle of serum, avoiding the loss of lyophilizate, and pipette exactly 5.0 ml of distilled / deionized water. Carefully close the bottle and dissolve the contents completely by occasional gentle swirling within 30 minutes.

Avoid the formation of foam.

Important: with the exception of **alkaline phosphatase**, all enzymes can be measured immediately. To reactivate the alkaline phosphatase, allow the reconstituted control to stand for one hour at +25°C.

Storage and stability

Store at 2-8°C. Stability of lyophilized control serum at 2-8°C: up to the stated expiration date.

Stability of the components in the reconstituted serum:

at +25°C	12 hours
at +4°C	5 days
at -20°C	1 month (when frozen once)

Stability of bilirubin in the reconstituted serum:

(when stored protected from light):

at +25°C	2 hours
at +4°C	6 hours
at -20°C	2 weeks (when frozen once)

The possible appearance of a slight green coloration has no effect on the recovery of the values. Store controls tightly capped when not in use.

Test procedure

Use supplied materials as provided.

Additional materials required: automatic micropipette, distilled/deionized water, general laboratory equipment.

Assay

Dispense the required volume into a sample cup and analyze in the same way as for the samples. The control should be run daily in parallel with the patient samples and after every calibration. The control intervals should be adapted to each laboratory's individual requirements. Each laboratory should establish QC procedures that conform with local, state and/or federal regulations or accreditation required.

Assigned values and ranges

Target value determination procedure: single assays in 7 external laboratories, with ≥ 3 independent series, under the direction of manufacturer department for external value assignment; this department also performed the evaluation. The declared value is the mean of all values obtained. The corresponding control range is calculated as the target value ± 3 SD, with the standard deviation (SD) being the value obtained from several target value determinations. Determinations were made using the testing reagents and analysis systems available at the time. Results must be within the defined ranges. Each laboratory should establish guidelines for corrective measures to be taken if values fall outside the range.

References

- 1) Occupational Safety and Health Standards: bloodborne pathogens. (29 CFR 1910.1030). Federal Register. July 1, 1998; 6:267-280
- 2) Council Directive (90/679/EEC). Official Journal of the European Communities. No. L374 from Dec. 31, 1990: 1-12.
- 3) Data on file

LTA srl
via Milano 15/F - 20060
Bussero (MI) - ITALY
Tel. ++39 02 95409034
Fax. ++39 02 95334185

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