

# INORGANIC PHOSPHORUS

## Colorimetric determination of inorganic phosphorus in biological liquids

### TEST SUMMARY

Phosphorus reacts with Ammonium molybdate forming a compound that absorbs ultra-violet light.

### SAMPLES

Serum, plasma immediately separated from cells.  
Diluted urine 1:25 with distilled water.  
Stability: 7 days at 4°C.

### REAGENTS

Sole reagent Ammonium Molybdate 0.8 M;  
Sulphuric Acid 0.6 M; surfactants,  
stabilizers.

Standard Inorganic phosphorus 4 mg/dl;  
stabilizers and conservatives.

### MATERIAL REQUIRED BUT NOT SUPPLIED

Normal laboratory equipment. Spectrophotometer UV/VIS with thermostatisation. Automatic Micropipette. Cuvette in optical glass or monouse in optical polystyrene. Physiologic solution.

### PRECAUTIONS

Reagent may contain not reactive and conservative components. It is opportune to avoid contacts with the skin and do not swallow.  
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

### REAGENTS PREPARATION

Reagents are ready to use, stable until expiration date on label, stored at 4-30°C.

### PROCEDURE

Kind of analysis: End Point  
Reading time: 2 minutes  
Wavelength: 340 nm (334-365)  
Temperature: R.T.  
Lightpath: 1 cm  
Zero: Blank Reagent  
Colour stability: 120 minutes

Reagents	Blank	Standard	Sample
Distilled water	10 µl	--	--
Standard	--	10 µl	--
Sample	--	--	10 µl
Sole Reagent	1000 µl	1000 µl	1000 µl

### CALCULATION

#### SERUM/PLASMA

Phosphorus (mg/dl) (A sample/A standard) x 4

Phosphorus (mM) (A sample /A standard) x 1.28

Phosphorus (mEq/l) (A sample /A standard) x 2.32

#### URINE (g/24h)

(A sample /A standard) x l of urine

### EXPECTED VALUES

	ADULTS	CHILDREN
Serum Plasma	2.8 - 4.6 mg/dl	4 - 6.5 mg/dl
	0.9 - 1.5 mM	1.3 - 2.1 mM
	1.6 - 2.7 mEq/l	2.3 - 3.8 mEq/l
Urine	0.5 - 1.1 g/24h	

Every laboratory should establish own reference intervals in accordance with own population.

### NOTES

- If the results are incompatible with clinical presentation, they have to be evaluated within a total clinical study.
- Only for IVD use.

### CALIBRATION/QUALITY CONTROL

It is suggested to perform an internal quality control. For this purpose the following control sera on human base are available on request:

**QN 0050 CH** 10 x 5 ml  
Control Sera normal values

**QP 0050 CH** 10 x 5 ml  
Control Sera pathological values

### TEST PERFORMANCE

#### Precision

Intra-assay (n = 20)	Mean (mg/dl)	SD (mg/dl)	CV%
Sample 1	3.2	0.032	1.01
Sample 2	4.3	0.032	0.75

Inter-assay (n = 20)	Mean (mg/dl)	SD (mg/dl)	CV%
Sample 1	3.21	0.044	1.39
Sample 2	4.305	0.039	0.92

#### Linearity

The method is linear up to 20 mg/dl.

#### Methods comparison

A comparison with a commercial available product gave the following results in a comparison on 50 samples:

Phosphorus LTA = x  
Phosphorus competitor = y  
n = 50

$$y = 0,99453x + 0,01144 \quad r = 0,9978$$

#### WASTE DISPOSAL

Product is intended for professional laboratories. Waste products must be handled as per relevant security cards and local regulations.

#### PACKAGING

**CODE CC01700 (600 TESTS)**  
Sole Reagent 6 x 100 ml (liquid)  
Standard 1 x 5 ml (liquid)









#### REFERENCES

Amodard – Clin. Chem. 18; 601 (1972).

#### MANUFACTURER

LTA s.r.l.  
Via Milano 15/F  
20060 Bussero (Milan) ITALY  
Tel: ++39 02 95409034  
Fax: ++39 02 95334185  
e-mail: info@ltaonline.it  
Website: http://www.ltaonline.it

### SYMBOLS

-  Only for IVD use
-  Lot of manufacturing
-  Code number
-  Storage temperature interval
-  Expiration date
-  Warning, read enclosed documents
-  Read the directions
-  Biological risk

Mod. 01.06 (ver. 4.3 - 08/07/2016)

