

XILOSE - CONTROL SET

Protein solutions to be used as Controls un Xilose determination in serum

INTENDED USE

Concentrations and activities of components are studied for an excellent use even in automatic system of clinical chemistry, even in manual procedures.

PRESERVATION AND STABILITY

Preserve the Controls at 2-8°C, stable until expiration date on the package.
Once opened controls are stable for 15 days at 2-8°C.

MATERIALS REQUIRED BUT NOT SUPPLIED

Automatic Micropipette, normal laboratory equipment.

REAGENTS PREPARATIONS

The controls are in liquid form and ready to use.
Bring controls at room temperature before use.

PRECAUTIONS

- Product contains Sodium Azide 0.095% as preservative.
- Avoid contact with skin and mucous membrane.
- Use the normal precautions for lab behaviour.
- For in vitro use only.

WASTE DISPOSAL

Product is intended for professional laboratories. Waste products must be handled as per relevant security cards and local regulations.

PACKAGING

CODE CC02630



Control Level 1 2 x 1 ml (liquid)
Control Level 2 2 x 1 ml (liquid)

REFERENCES

Archives data.

CALIBRATORS VALUES









LOT C165.18*  2020/06

CALIBRATORS VALUES		CONTROL LEVEL 1		CONTROL LEVEL 2	
		LOT G157.18*  2020/06		LOT H157.18*  2020/06	
METHOD	U	VALUE	RANGE	VALUE	RANGE
Enzymatic U.V.	mg/dl	15,23	12,18 – 18,28	51,53	41,22 – 61,84

MANUFACTURER

LTA s.r.l.
Via Milano 15/F
20060 Bussero (Milan) ITALY
Tel: +39 02 95409034
Fax: +39 02 95334185
e-mail: info@LTAonline.it
Website: www.LTAonline.it

SYMBOLS

-  Only for IVD use
-  Lot of manufacturing
-  Code number
-  Storage temperature interval
-  Expiration date (year – month)
-  Warning, read enclosed documents
-  Read the directions
-  Biological risk

Mod. 01.06 (ver. 2.0 – 07/06/2018)

