

APOLIPOPROTEIN A₁/B CALIBRATOR (CALIBRATOR FOR IMMUNOCHEMICAL REACTIONS)

Titred human serum useful as calibrator or control in immuno turbidimetric and radial immunodiffusion reactions

PRINCIPLE

The evaluation of suitable material in the sequence analysis, allows to verify the global inaccuracy of procedure adopted in. The use of human serum as a control material, guarantees commutability of results obtained in the Quality Control process with similar samples.

REAGENTS

Treated human serum ready to use.

STORAGE AND STABILITY

The vial is stored at 2-8°C and is stable until expiration date on label if not contaminated during use.

Before opening wait and bring the vial at room temperature, and shake delicately avoiding foam formation.

MATERIALS REQUIRED BUT NOT PROVIDED

Automatic micropipette, deionized or distilled water, general laboratory equipment.

CALIBRATION CURVE

To obtain a calibration curve, dilute the calibrator with saline solution at different concentrations.

Calibrators must be treated as if they were samples and test as indicated by the analytical procedure required by the reagent used.

REFERENCE INTERVALS

The reference intervals are defined taking account of preparation variability, reconstitution products and equipment, by means of a statistical process of the data collected.

The component values are reported in the attached table.

PRECAUTIONS

- The blood of donors of biological substances used for products preparation, has been tested and found negative for HbsAg and antibodies anti-HIV and HCV.
- Considering that any known test can guarantee with absolute certainty that products derived from human blood, don't have risk of Hepatitis or HIV infections, the material should be handled just as carefully as a patient sample.

NOTE

- Inaccurate results, changes of color reagent or loss of vacuum in the vials could indicate a deterioration of the product. In any case, poor performance of the calibrator may also be due to other factors related to the dosage system.
- For in vitro diagnostic use only.

WASTE DISPOSAL

Product is intended for professional laboratories.

Waste products must be handled as per relevant security cards and local regulations.

PACKAGING









CODE IC00300

Apolipoprotein A₁/B Calibrator 1 x 1 ml (liquid)

MANUFACTURER

LTA s.r.l.
Via Milano 15/F
20060 Bussero (Milano) - ITALY
tel. ++39 02 95409034
fax. ++39 02 95334185
e-mail. info@ltaonline.it
website. <http://www.ltaonline.it>

SYMBOLS

-  Only for IVD use
-  Lot of manufacturing
-  Code number
-  Storage temperature interval
-  Expiration date (year, month)
-  Warning, read enclosed documents
-  Read the directions
-  Biological risk

Mod. 01.06 (ver. 1.0 – 05/12/2009)

