

# KAPPA / LAMBDA TOTAL LIGHT CHAINS CALIBRATOR (CALIBRATOR FOR IMMUNOCHEMICAL REACTIONS)

**Titred human serum useful as calibrator or control in immuno turbidimetric and radial immunodiffusion reactions for urine tests**

## PRINCIPLE

The evaluation of suitable material in the sequence analysis, allows to verify the global inaccuracy of procedure adopted in.

The use of human serum as a control material, guarantees commutability of results obtained in the Quality Control process with similar samples.

The serum concentrations are specially designed to be optimal for use in testing of the urine, and are calibrated against ERM-DA470k/IFCC.

(SEE POINT VALUES CALIBRATION).

## REAGENTS

Treated and titred human serum ready to use.

## STORAGE AND STABILITY

Calibrator is stored at 2-8°C and is stable until expiration date on label if not contaminated during use.

Before opening wait and bring the vial at room temperature, and shake delicately avoiding foam formation.

## VALUES CALIBRATION

The calibration values are specified for each lot and are indicated in the attached sheet values.

Kappa and Lambda values are calculated using the concentrations of individual proteins to IgG, IgA, and IgM determined by the average of the tests performed with immunoturbidimetric and radial immuno-diffusion methods using as a calibrator ERM-DA470k/IFCC and applying the following formula published by M.M. Lievens in J Clin Chem Clin Biochem 1989;27:519-23:

$$\text{Kappa} = [\text{IgG}] \times 0,1983 + [\text{IgA}] \times 0,1710 + [\text{IgM}] \times 0,0975$$
$$\text{Lambda} = [\text{IgG}] \times 0,1054 + [\text{IgA}] \times 0,1206 + [\text{IgM}] \times 0,0305$$

The extended uncertainty is calculated with  $k = 2$  corresponds to a confidence level of 95%.

## CALIBRATION CURVE

To obtain a calibration curve, dilute the calibrator with saline solution at different concentrations.

Calibrators must be treated as if they were samples and test as indicated by the analytical procedure required by the reagent used.

## REFERENCE INTERVALS

Reference ranges for protein for healthy individuals.

### *In Serum*

Kappa total light chains	140 – 380 mg/dl
Lambda total light chains	90 – 245 mg/dl

### *In Urine*

Kappa total light chains	< 4 mg/dl
Lambda total light chains	< 4 mg/dl

Each laboratory should establish appropriate reference intervals related to its population.

## PRECAUTIONS

- Contains Sodium Azide 0.09 g/l.
- Each serum or plasma donor unit used in this product has been tested and found to be negative for HbsAg, HCV, HIV antibody, HCV-PCR and HIV-PCR by an FDA approved method.
- Considering that no known test can offer complete certainty that products derived from human blood will not transmit hepatitis, HIV or other infectious agents, must be taken all the

necessary precautions for handling potentially infectious materials.

## NOTE

- Inaccurate results, or changes of color reagent could indicate a deterioration of the product. In any case, poor performance of the calibrator may also be due to other factors related to the dosage system.
- For in vitro diagnostic use only.

## WASTE DISPOSAL

Waste products must be handled as per relevant security cards and local regulations.

## PACKAGING

### CODE IC00900

Calibrator 1 x 1 ml (liquid)


## MANUFACTURER


LTA s.r.l.  
Via Milano 15/F  
20060 Bussero (Milano) - ITALY  
tel. ++39 02 95409034  
fax. ++39 02 95334185  
e-mail. info@ltaonline.it  
website. http://www.ltaonline.it

## REFERENCES

1. Blirup-Jensen S, Johnson AM, Larsen M. Protein standardization V: value transfer. A practical protocol for the assignment of serum protein values from a reference material to a Target Material. Clin Chem Lab Med 2008;46:1470–79.
2. Zegers I, Keller T, Schreiber W, Sheldon J, Albertini R, Blirup-Jensen S, Johnson M, Trapmann S, Emons H, Merlini G, Schimmel H. Characterization of the New Serum Protein Reference Material ERM-DA470k/IFCC: Value Assignment by Immunoassay. Clin Chem 2010;56:12:1880–88.
3. Blirup-Jensen S, Johnson AM, Larsen M. Protein standardization IV: value transfer. Procedure for the assignment of serum protein values from a reference preparation to a target material. Clin Chem Lab Med 2001;39:1110–22.
4. M.M. Lievens in J Clin Chem Clin Biochem 1989;27:519-23


## SYMBOLS


 Only for IVD use


 Lot of manufacturing

 Code number

 Storage temperature interval

 Expiration date (year, month)

 Warning, read enclosed documents

 Read the directions

 Biological risk

Mod. 01.06 (ver. 1.0 – 15/09/2016)

